

### **Test Definition: CD279**

CD279 (PD-1) Immunostain, Technical Component Only

# **Reporting Title:** CD279 (PD-1) IHC, Tech Only **Performing Location:** Rochester

#### **Ordering Guidance:**

This test includes only technical performance of the stain (no pathologist interpretation is performed). If diagnostic consultation by a pathologist is required order PATHC / Pathology Consultation.

#### **Shipping Instructions:**

Attach the green pathology address label and the pink Immunostain Technical Only label included in the kit to the outside of the transport container.

#### **Specimen Requirements:**

Supplies: Immunostain Technical Only Envelope (T693)
Specimen Type: Tissue
Container/Tube: Immunostain Technical Only Envelope
Preferred: 2 unstained positively charged glass slide (25- x 75- x 1-mm) per test ordered; sections 4-microns thick.
Acceptable: Formalin-fixed, paraffin-embedded (FFPE) tissue block

#### Forms:

If not ordering electronically, complete, print, and send a <u>Immunohistochemical (IHC)/In Situ Hybridization (ISH) Stains</u> <u>Request</u> (T763) with the specimen.

Specimen Type	Temperature	Time	Special Container
TECHONLY Ambient (preferred)			
	Refrigerated		

#### **Result Codes:**

Result ID	Reporting Name	Туре	Unit	LOINC®
70730	CD279 (PD-1) IHC, Tech Only	Alphanumeric		Bill only; no result

LOINC<sup>®</sup> and CPT codes are provided by the performing laboratory.

#### Supplemental Report:

No

#### **CPT Code Information:**

88342-TC, primary 88341-TC, if additional IHC

#### **Reflex Tests:**

Т	est Id	Reporting Name	CPT Units	CPT Code	Always Performed	Available Separately



## **Test Definition: CD279**

CD279 (PD-1) Immunostain, Technical Component Only

IHTOI	IHC Initial, Tech Only	1	88342	No	No
IHTOA	IHC Additional, Tech Only	1	88341	No	No

**Reference Values:**