

Reporting Title: FUS IHC, Tech Only
Performing Location: Rochester

Ordering Guidance:
This test includes only technical performance of the stain (no pathologist interpretation is performed). If diagnostic consultation by a pathologist is required order PATHC / Pathology Consultation.

Shipping Instructions:
Attach the green "Attention Pathology" address label (T498) and the pink Immunostain Technical Only label included in the kit to the outside of the transport container.

Specimen Requirements:
Supplies: Immunostain Technical Only Envelope (T693)
Specimen Type: Tissue
Container/Tube: Immunostain Technical Only Envelope
Preferred: 2 Unstained positively charged glass slide (25- x 75- x 1-mm) per test ordered; sections 4-microns thick
Acceptable: Formalin-fixed, paraffin-embedded tissue block

Forms:
If not ordering electronically, complete, print, and send a [Immunohistochemical \(IHC\)/In Situ Hybridization \(ISH\) Stains Request](#) (T763) with the specimen.

| Specimen Type | Temperature | Time | Special Container |
|---------------|---------------------|------|-------------------|
| TECHONLY | Ambient (preferred) | | |
| | Refrigerated | | |

Result Codes:

| Result ID | Reporting Name | Type | Unit | LOINC® |
|-----------|--------------------|--------------|------|----------------------|
| 70747 | FUS IHC, Tech Only | Alphanumeric | | Bill only; no result |

LOINC® and CPT codes are provided by the performing laboratory.

Supplemental Report:
No

CPT Code Information:
88342-TC, primary
88341-TC, if additional IHC

Reflex Tests:

| Test Id | Reporting Name | CPT Units | CPT Code | Always Performed | Available Separately |
|---------|------------------------|-----------|----------|------------------|----------------------|
| IHTOI | IHC Initial, Tech Only | 1 | 88342 | No | No |

| | | | | | |
|-------|---------------------------|---|-------|----|----|
| IHTOA | IHC Additional, Tech Only | 1 | 88341 | No | No |
|-------|---------------------------|---|-------|----|----|