

Test Definition: KIMEL

Ki67 + Melan A Immunostain, Technical Component Only

Reporting Title: Ki67+Melan A IHC, Tech Only **Performing Location:** Rochester

Ordering Guidance:

This test includes only technical performance of the stain (no pathologist interpretation is performed). If diagnostic consultation by a pathologist is required order PATHC / Pathology Consultation.

Shipping Instructions:

Attach the green "Attention Pathology" address label (T498) and the pink Immunostain Technical Only label included in the kit to the outside of the transport container.

Specimen Requirements:

Supplies: Immunostain Technical Only Envelope (T693)
Specimen Type: Tissue
Container/Tube: Immunostains Technical Only Envelope
Preferred: 2 Unstained positively charged glass slides (25- x 75- x 1-mm) per test ordered; sections 4-microns thick.
Acceptable: Formalin-fixed, paraffin-embedded tissue block

Forms:

If not ordering electronically, complete, print, and send a <u>Immunohistochemical (IHC)/In Situ Hybridization (ISH) Stains</u> <u>Request</u> (T763) with the specimen.

| Specimen Type | Temperature | Time | Special Container |
|---------------|---------------------|------|-------------------|
| TECHONLY | Ambient (preferred) | | |
| | Refrigerated | | |

Result Codes:

| Result ID | Reporting Name | Туре | Unit | LOINC® |
|-----------|-----------------------------|--------------|------|----------------------|
| 70795 | Ki67+Melan A IHC, Tech Only | Alphanumeric | | Bill only; no result |

LOINC® and CPT codes are provided by the performing laboratory.

Supplemental Report:

No

CPT Code Information:

88344-TC

Reflex Tests:

| Test | t Id | Reporting Name | CPT Units | CPT Code | Always Performed | Available Separately |
|------|------|--------------------------|-----------|----------|------------------|----------------------|
| IHM | ITO | IHC Multiplex, Tech Only | 1 | 88344 | No | No |



Test Definition: KIMEL

Ki67 + Melan A Immunostain, Technical Component Only