

Test Definition: P40NA

p40 + Napsin A Immunostain, Technical Component Only

Reporting Title: p40 + Napsin A IHC, Tech Only **Performing Location:** Rochester

Ordering Guidance:

This test includes only technical performance of the stain (no pathologist interpretation is performed). If diagnostic consultation by a pathologist is required order PATHC / Pathology Consultation.

Shipping Instructions:

Attach the green "Attention Pathology" address label (T498) and the pink Immunostain Technical Only label included in the kit to the outside of the transport container.

Specimen Requirements:

Supplies: Immunostain Technical Only Envelope (T693)
Specimen Type: Tissue
Container/Tube: Immunostains Technical Only Envelope
Preferred: 2 Unstained positively charged glass slide (25- x 75- x 1-mm) per test ordered; sections 4-microns thick.
Acceptable: Formalin-fixed, paraffin-embedded tissue block

Forms:

If not ordering electronically, complete, print, and send a <u>Immunohistochemical (IHC)/In Situ Hybridization (ISH) Stains</u> <u>Request</u> (T763) with the specimen.

Specimen Type	Temperature	Time	Special Container
TECHONLY Ambient (preferred)			
	Refrigerated		

Result Codes:

Result ID	Reporting Name	Туре	Unit	LOINC®
70838	p40 + Napsin A IHC, Tech Only	Alphanumeric		Bill only; no result

LOINC[®] and CPT codes are provided by the performing laboratory.

Supplemental Report:

No

CPT Code Information:

88344-TC

Reflex Tests:

Test	t Id	Reporting Name	CPT Units	CPT Code	Always Performed	Available Separately
IHM	ITO	IHC Multiplex, Tech Only	1	88344	No	No



Test Definition: P40NA

p40 + Napsin A Immunostain, Technical Component Only