

Test Definition: PLAP

Placental Alkaline Phosphatase (PLAP) Immunostain, Technical Component Only

Reporting Title: PLAP IHC, Tech Only

Performing Location: Rochester

Ordering Guidance:

This test includes only technical performance of the stain (no pathologist interpretation is performed). If diagnostic consultation by a pathologist is required order PATHC / Pathology Consultation.

Shipping Instructions:

Attach the green "Attention Pathology" address label (T498) and the pink Immunostain Technical Only label included in the kit to the outside of the transport container.

Specimen Requirements:

Supplies: Immunostain Technical Only Envelope (T693)

Specimen Type: Tissue

Container/Tube: Immunostains Technical Only Envelope

Preferred: 2 Unstained positively charged glass slide (25- x 75- x 1-mm) per test ordered; sections 4-microns thick.

Acceptable: Formalin-fixed, paraffin-embedded tissue block

Forms:

If not ordering electronically, complete, print, and send a <u>Immunohistochemical (IHC)/In Situ Hybridization (ISH) Stains</u>
Request (T763) with the specimen.

Specimen Type	Temperature	Time	Special Container
TECHONLY	Ambient (preferred)		
	Refrigerated		

Result Codes:

Result ID	Reporting Name	Туре	Unit	LOINC®
70851	PLAP IHC, Tech Only	Alphanumeric		Bill only; no result

LOINC® and CPT codes are provided by the performing laboratory.

Supplemental Report:

No

CPT Code Information:

88342-TC, primary 88341-TC, if additional IHC

Reflex Tests:

Test Id	Reporting Name	CPT Units	CPT Code	Always Performed	Available Separately
---------	----------------	-----------	----------	------------------	----------------------



Test Definition: PLAP

Placental Alkaline Phosphatase (PLAP) Immunostain, Technical Component Only

IHTOI	IHC Initial, Tech Only	1	88342	No	No
IHTOA	IHC Additional, Tech Only	1	88341	No	No