

Reporting Title: Kappa/Lambda ISH, Tech Only
Performing Location: Rochester

Ordering Guidance:
This test includes only technical performance of the stain (no pathologist interpretation is performed). If diagnostic consultation by a pathologist is required order PATHC / Pathology Consultation.

Shipping Instructions:
Attach the green "Attention Pathology" address label (T498) and the pink Immunostain Technical Only label included in the kit to the outside of the transport container.

Specimen Requirements:
Supplies: Immunostain Technical Only Envelope (T693)
Specimen Type: Tissue
Container/Tube: Immunostain Technical Only Envelope
Preferred: 5 Unstained positively charged glass slides (25- x 75- x 1-mm) per test ordered; sections 4-microns thick.
Acceptable: Formalin-fixed, paraffin-embedded tissue block

Forms:
If not ordering electronically, complete, print, and send a [Immunohistochemical \(IHC\)/In Situ Hybridization \(ISH\) Stains Request](#) (T763) with the specimen.

Specimen Type	Temperature	Time	Special Container
TECHONLY	Ambient (preferred)		
	Refrigerated		

Result Codes:

Result ID	Reporting Name	Type	Unit	LOINC®
70792	Kappa/Lambda ISH, Tech Only	Alphanumeric		Bill only; no result

LOINC® and CPT codes are provided by the performing laboratory.

Supplemental Report:
No

CPT Code Information:
88365-TC, primary
88364-TC, if additional ISH

Reflex Tests:

Test Id	Reporting Name	CPT Units	CPT Code	Always Performed	Available Separately
---------	----------------	-----------	----------	------------------	----------------------

Test Definition: KLISH

Kappa and Lambda Light Chain mRNA, In Situ
Hybridization (ISH) Technical Component Only

ISTOI	ISH Initial, Tech Only	1	88365	No	No
ISTOA	ISH Additional, Tech Only	1	88364	No	No