

# **Test Definition: FALPX**

Alprazolam (Xanax)

Reporting Title: Alprazolam (Xanax)

Performing Location: Medtox Laboratories, Inc.

## **Specimen Requirements:**

Submit only 1 of the following specimens:

#### Serum

Draw blood in a plain, red-top tube(s), serum gel tube is not acceptable. Spin down and send 2 mL of serum refrigerated in a plastic vial.

#### Plasma

Draw blood in a green-top sodium heparin tube(s), plasma gel tube is not acceptable. Spin down and send 2 mL of plasma refrigerated in a plastic vial.

Specimen Type	Temperature	Time	Special Container
Varies	Refrigerated (preferred)	14 days	
	Frozen	180 days	
	Ambient	72 hours	

### **Result Codes:**

Ī	Result ID	Reporting Name	Туре	Unit	LOINC®
	FALPX	Alprazolam (Xanax)	Alphanumeric		59611-4

LOINC® and CPT codes are provided by the performing laboratory.

# **Supplemental Report:**

No

#### **CPT Code Information:**

80346

### **Reference Values:**

5 - 25 ng/mL

Reporting Limit: 2.0 ng/mL