

Acetoacetate, Serum

Reporting Title: Acetoacetate, Serum **Performing Location:** NMS Labs

Specimen Requirements:

Must submit one specimen per order. Specimens cannot be shared between multiple orders.

Specimen Type: Serum

Collection Container/Tube: Red

Submission Container/Tube: Plastic vial

Specimen Volume: 3 mL

Collection Instructions:

1. Draw blood in a plain, red-top tube(s). Serum gel tube is not acceptable.

2. Centrifuge and send 3 mL of serum frozen in a plastic, preservative-free vial.

Specimen Type	Temperature	Time	Special Container
Serum Red	Frozen (preferred)	4 days	
	Refrigerated	4 days	

Result Codes:

Result ID	Reporting Name	Туре	Unit	LOINC®
Z3524	Acetoacetate	Alphanumeric		Not Provided

LOINC[®] and CPT codes are provided by the performing laboratory.

Supplemental Report:

No

CPT Code Information:

82010

Reference Values:

Reporting limit determined each analysis.

Acetoacetate: None Detected mcg/mL -Normal range for adults: 5-30 mcg/mL