

Sinemet, Serum

**Reporting Title:** Sinemet, Serum **Performing Location:** NMS Labs

# **Specimen Requirements:**

\*\*\*Must submit one specimen per order. Specimens cannot be shared between multiple orders.\*\*\*

Specimen Type: Serum

Collection Container/Tube: Red top

Submission Container/Tube: Plastic vial

Specimen Volume: 1 mL

#### **Collection Instructions:**

1. Draw blood in a plain, red-top tube(s). Serum gel tube is not acceptable.

2. Centrifuge and send 1 mL of serum frozen in a preservative-free plastic vial. **Specimen must be frozen immediately**, or results will be compromised.

| Specimen Type | Temperature | Time   | Special Container |
|---------------|-------------|--------|-------------------|
| Serum Red     | Frozen      | 6 days |                   |

## **Result Codes:**

| Result ID | Reporting Name | Туре         | Unit | LOINC® |
|-----------|----------------|--------------|------|--------|
| Z3535     | Levodopa       | Alphanumeric |      | 9385-6 |
| Z3536     | Carbidopa      | Alphanumeric |      | 9328-6 |

LOINC<sup>®</sup> and CPT codes are provided by the performing laboratory.

# Supplemental Report:

No

## **CPT Code Information:**

80299

## **Reference Values:**

Reporting limit determined each analysis.

Levodopa Units: mcg/mL

Carbidopa Units: mcg/mL