

# **Test Definition: FCLMS**

Clonazepam and Metabolite, Serum

## **Reporting Title:** Clonazepam and Metabolite, S **Performing Location:** NMS Labs

### **Specimen Requirements:**

\*\*\*Must submit one specimen per order. Specimens cannot be shared between multiple orders.\*\*\*

Specimen Type: Serum

Collection Container/Tube: Red top

Submission Container/Tube: Plastic, preservative-free vial

Specimen Volume: 2 mL

#### **Collection Instructions:**

- 1. Draw blood in a red top tube(s). Serum gel tube is not acceptable.
- 2. Centrifuge and send 2 mL of serum refrigerated in a plastic, preservative-free vial.

Specimen Type	Temperature	Time	Special Container
Serum Red	Refrigerated (preferred)	14 days	
	Frozen	30 days	

#### Result Codes:

Result ID	Reporting Name	Туре	Unit	LOINC®
Z6336	Clonazepam	Alphanumeric		3494-2
Z6337	7-Amino Clonazepam	Alphanumeric		28059-4

LOINC<sup>®</sup> and CPT codes are provided by the performing laboratory.

#### Supplemental Report:

No

#### **CPT Code Information:**

80346 G0480 (if appropriate)

#### **Reference Values:**

Reporting limit(s) determined each analysis. Clonazepam: None Detected ng/mL 7-Amino Clonazepam: None Detected ng/mL