

Reporting Title: Meconium Methadone Scr/Confirm
Performing Location: Medtox Laboratories, Inc.

Specimen Requirements:

Container/Tube: Sterile collection container

Specimen Volume: 3 g

Collection Instructions: Collect specimen into the same sterile collection container until 3 g of meconium have been collected or until the first milk stool appears. When at least 3 g of meconium has been collected, tightly screw on the cap of the collection vial and send specimen ambient.

Note: Specimens from different voids may be pooled if necessary.

| Specimen Type | Temperature | Time | Special Container |
|---------------|---------------------|----------|-------------------|
| Meconium | Ambient (preferred) | 14 days | |
| | Frozen | 365 days | |
| | Refrigerated | 14 days | |

Result Codes:

| Result ID | Reporting Name | Type | Unit | LOINC® |
|-----------|----------------|--------------|------|---------|
| Z5603 | Methadone | Alphanumeric | | 43595-8 |

LOINC® and CPT codes are provided by the performing laboratory.

Supplemental Report:

No

CPT Code Information:

80307

Reflex Tests:

| Test Id | Reporting Name | CPT Units | CPT Code | Always Performed | Available Separately |
|---------|----------------------------|-----------|----------|------------------|----------------------|
| FMTHM | Meconium Methadone Confirm | 1 | 80358 | No | No |

Result Codes for Reflex Tests:

| Test ID | Result ID | Reporting Name | Type | Unit | LOINC® |
|---------|-----------|-----------------------------|--------------|------|---------|
| FMTHM | Z5604 | Methadone | Alphanumeric | | 29159-1 |
| FMTHM | Z5605 | Methadone Metabolite (EDDP) | Alphanumeric | | 73562-1 |

Reference Values:

Initial presumptive testing by immunoassay at the following testing threshold:

Methadone: 50 ng/gm

Presumptive positive confirmed by definitive chromatography mass spectrometry (GC/MC or LC/MS-MS)