

Reporting Title: ANA, Synovial Fluid

Performing Location: Esoterix Endocrinology

Ordering Guidance:

Specimen Requirements:

Specimen Type: Synovial fluid

Collection Container/Tube: Sterile container

Submission Container/Tube: Plastic vial

Specimen Volume: 1 mL

Collection Instructions:

- 1. Collect synovial fluid and transfer into a plastic vial and ship refrigerate.
- 2. Indicate the specimen source.

| Specimen Type | Temperature              | Time    | Special Container |
|---------------|--------------------------|---------|-------------------|
| Varies        | Refrigerated (preferred) | 14 days |                   |
|               | Frozen                   | 60 days |                   |
|               | Ambient                  | 7 days  |                   |

Ask at Order Entry (AOE) Questions:

| Test ID | Question ID | Description       | Type       | Reportable |
|---------|-------------|-------------------|------------|------------|
| FFANA   | Z5710       | ANA Specimen Type | Plain Text | Yes        |

Result Codes:

| Result ID | Reporting Name      | Type         | Unit | LOINC®       |
|-----------|---------------------|--------------|------|--------------|
| Z5710     | ANA Specimen Type   | Alphanumeric |      | Not Provided |
| Z5711     | ANA, Synovial Fluid | Alphanumeric |      | In Process   |

LOINC® and CPT codes are provided by the performing laboratory.

Supplemental Report:

No

CPT Code Information:

86038

86039 (if appropriate)

Reflex Tests:

| Test Id | Reporting Name                       | CPT Units | CPT Code | Always Performed | Available Separately |
|---------|--------------------------------------|-----------|----------|------------------|----------------------|
| FPANA   | ANA Titer and Pattern Synovial Fluid |           |          | No               | No                   |

Result Codes for Reflex Tests:

| Test ID | Result ID | Reporting Name            | Type         | Unit | LOINC®  |
|---------|-----------|---------------------------|--------------|------|---------|
| FPANA   | Z5712     | Note:                     | Alphanumeric |      | 13067-4 |
| FPANA   | Z5729     | Homogeneous Pattern       | Alphanumeric |      | 13067-4 |
| FPANA   | Z5730     | Nucleolar Pattern         | Alphanumeric |      | 13067-4 |
| FPANA   | Z5731     | Speckled Pattern          | Alphanumeric |      | 13067-4 |
| FPANA   | Z5732     | Centromere Pattern        | Alphanumeric |      | 13067-4 |
| FPANA   | Z5733     | Spindle Apparatus Pattern | Alphanumeric |      | 13067-4 |
| FPANA   | Z5734     | Nuclear Membrane Pattern  | Alphanumeric |      | 13067-4 |
| FPANA   | Z5735     | Midbody Pattern           | Alphanumeric |      | 13067-4 |
| FPANA   | Z5736     | Nuclear Dot Pattern       | Alphanumeric |      | 13067-4 |
| FPANA   | Z5737     | PCNA Pattern              | Alphanumeric |      | 13067-4 |
| FPANA   | Z5738     | Centriole Pattern         | Alphanumeric |      | 13067-4 |

Reference Values:

ANA Titer: <1:10  
ANA Pattern: No Pattern