

Reporting Title: Rifampin Level
Performing Location: National Jewish Health

Specimen Requirements:

Container/Tube: Red Top
Preferred: Red top tube
Acceptable: Serum gel tube
Specimen Volume: 2 mL

Collection Instructions: Draw blood in a plain, red-top tube(s). Separate serum from cells immediately by centrifugation and aliquot into a polypropylene or similar plastic tube. Send 2 mL of serum frozen in plastic vial.

- Note:
- The following information is required:
 - Specimen Type (source)
 - Dose (specify PO, IV, IM)
 - Date and time of last dose (for IV start/end time)
 - If the time of last dose and the blood draw are not accurately recorded, accurate interpretation of the concentration is not possible.

Specimen Type	Temperature	Time	Special Container
Serum	Frozen	365 days	

Ask at Order Entry (AOE) Questions:

Test ID	Question ID	Description	Type	Reportable
FRIFA	Z5803	Dose	Plain Text	Yes
FRIFA	Z5804	Date and time of last dose	Plain Text	Yes
FRIFA	Z5865	Specimen Type:	Plain Text	Yes

Result Codes:

Result ID	Reporting Name	Type	Unit	LOINC®
Z5803	Dose	Alphanumeric	mL	Not Provided
Z5804	Date and Time of Last Dose	Alphanumeric	mL	Not Provided
Z5778	Rifampin Level	Alphanumeric	mL	4021-2
Z5840	Comment:	Alphanumeric	mL	Not Provided
Z5865	Specimen Type:	Alphanumeric	mL	Not Provided

LOINC® and CPT codes are provided by the performing laboratory.

Supplemental Report:

No

CPT Code Information:
80299

Reference Values:
mcg/mL