

Test Definition: FAS1A

Anti-SAE1 Ab, IgG

Reporting Title: Anti-SAE1 Ab, IgG
Performing Location: Esoterix Endocrinology

Specimen Requirements:

Collection Container/Tube: 5 mL Red **Submission Container/Tube:** Plastic vial

Specimen Volume: 1 mL

Acceptable: SST

Collection Instructions: Draw blood in a plain red-top tube(s), serum gel tube is acceptable. Spin down and send 1 mL of

serum refrigerated in a plastic vial.

Specimen Type	Temperature	Time	Special Container
Serum	Refrigerated (preferred)	14 days	
	Frozen	60 days	
	Ambient	7 days	

Result Codes:

Result ID	Reporting Name	Туре	Unit	LOINC®
FAS1A	Anti-SAE1 Ab, IgG	Alphanumeric		101220-2

LOINC® and CPT codes are provided by the performing laboratory.

Supplemental Report:

No

CPT Code Information:

83516

Reference Values:

Reference Range: <20

Interpretation:

Negative: <20 units

Weak Positive: 20-39 units Moderate Positive: 40-80 units Strong Positive: >80 units