

Test Definition: FFTAC

Thrombin-Antithrombin Complex

Reporting Title: Thrombin Antithrombin Complex

Performing Location: Esoterix Coagulation

Specimen Requirements: Specimen Type: Citrated plasma

Collection Container/Tube: Light-blue top (citrate)

Specimen Volume: 2 mL

Collection Instructions: Draw blood in a light blue-top (Sodium citrate) tube(s). Spin down and send 2 mL citrated

plasma frozen in a plastic vial.

Note: Separate specimens must be submitted when multiple tests are ordered.

Specimen Type	Temperature	Time	Special Container	
Plasma Na Cit	Frozen	180 days		

Result Codes:

Result ID	Reporting Name	Туре	Unit	LOINC®
FFTAC	Thrombin Antithrombin Complex	Alphanumeric		14182-0

LOINC® and CPT codes are provided by the performing laboratory.

Supplemental Report:

No

CPT Code Information:

83520

Reference Values:

<4.3 ng/mL

Pre-analytical conditions such as a difficult draw may spuriously increase test results.