

Reporting Title: Posaconazole (POS) Susceptibility
Performing Location: U.T. Health Science Ctr San Antonio

Specimen Requirements:

- Container/Tube:** Agar slant or other appropriate media
Specimen Volume: Organism in pure culture
Complete and submit with specimen:
- 1. Fungus Test Lab Request Form

Collection Instructions:

- 1. Organism must be in pure culture, actively growing.
- 2. Place specimen in a large infectious container (T146) and label as etiologic agent/infectious substance.

Specimen Type	Temperature	Time	Special Container
Varies	Ambient		

Result Codes:

Result ID	Reporting Name	Type	Unit	LOINC®
FPPOS	Posaconazole (POS) Susceptibility	Alphanumeric		Not Provided

LOINC® and CPT codes are provided by the performing laboratory.

Supplemental Report:

No

CPT Code Information:

- 87999 - Unlisted Microbiology Procedure (Refer to patient report to apply the appropriate CPT code in place of this unlisted CPT code based on specimen organism.)
87188 - mould - MIC microdilution or agar dilution (if appropriate)
87186 - yeast - MIC microdilution or agar dilution (if appropriate)