

# Test Definition: FLUNP

Influenza Virus Type A and Type B, and  
Respiratory Syncytial Virus (RSV), Molecular  
Detection, PCR, Nasopharyngeal Swab

**Reporting Title:** Influenza A/B and RSV, PCR, NP Swab  
**Performing Location:** Jacksonville

**Specimen Requirements:**

**Specimen Type:** Nasopharyngeal swab  
**Container/Tube:** Sterile container with viral transport media  
**Collection Instructions:**

- 1. Collect specimen by swabbing back and forth over mucosa surface to maximize recovery of cells.
- 2. Swab must be placed in viral transport media (for example, M4-RT, M4, or M4 media).

**Forms:**

Specimen Type	Temperature	Time	Special Container
Swab	Refrigerated (preferred)	7 days	
	Frozen	7 days	

**Result Codes:**

Result ID	Reporting Name	Type	Unit	LOINC®
SS019	Specimen Source	Alphanumeric		31208-2
35975	Influenza A, PCR	Alphanumeric		76078-5
35976	Influenza B, PCR	Alphanumeric		76080-1
35977	Respiratory Syncytial Virus, PCR	Alphanumeric		76089-2

LOINC® and CPT codes are provided by the performing laboratory.

**Supplemental Report:**

No

**CPT Code Information:**

87631

**Reference Values:**

Negative