

Test Definition: VIP

Vasoactive Intestinal Polypeptide, Plasma

Reporting Title: Vasoactive Intestinal Polypeptide,P

Performing Location: Rochester

Ordering Guidance:

This test should not be requested on patients who have recently received radioactive material.

Specimen Requirements:

Patient Preparation: Fasting (8 hours)

Collection Container/Tube: Lavender top (EDTA)

Submission Container/Tube: Plastic vial

Specimen Volume: 1 mL

Collection Instructions: Centrifuge and aliquot plasma into a plastic vial. Freeze immediately.

Forms:

If not ordering electronically, complete, print, and send an Oncology Test Request (T729) with the specimen.

Specimen Type	Temperature	Time	Special Container	
Plasma EDTA	Frozen	90 days		

Result Codes:

Result ID	Reporting Name	Туре	Unit	LOINC®
8150	Vasoactive Intestinal Polypeptide,P	Numeric	pg/mL	3125-2

LOINC® and CPT codes are provided by the performing laboratory.

Supplemental Report:

No

CPT Code Information:

84586

Reference Values:

<75 pg/mL