

Reporting Title: Methadone Confirmation, U

Performing Location: Rochester

Ordering Guidance:

- 1. For situations where chain of custody is required, a Chain-of-Custody Kit (T282) is available. For chain-of-custody testing, order MTDNX / Methadone Confirmation, Chain of Custody, Random, Urine.
- 2. Additional drug panels and specific requests are available. Call 800-533-1710 or 507-266-5700.
- 3. If urine creatinine is required or adulteration of the sample is suspected, order ADULT / Adulterants Survey, Random, Urine.

Specimen Requirements:

- Supplies: Urine Tubes, 10 mL (T068)
- Container/Tube: Plastic, 10-mL urine tube
- Specimen Volume: 10 mL
- Collection Instructions:

- 1. Collect a random urine specimen.
- 2. No preservative.

Forms:

If not ordering electronically, complete, print, and send a [Therapeutics Test Request](#) (T831) with the specimen.

Specimen Type	Temperature	Time	Special Container
Urine	Refrigerated (preferred)	28 days	
	Frozen	28 days	
	Ambient	14 days	

Result Codes:

Result ID	Reporting Name	Type	Unit	LOINC®
56028	EDDP-by GC-MS	Alphanumeric	ng/mL	58428-4
83129	Methadone-by GC-MS	Alphanumeric	ng/mL	16246-1
21107	Methadone Interpretation	Alphanumeric		69050-3
21110	Chain of Custody	Alphanumeric		77202-0

LOINC® and CPT codes are provided by the performing laboratory.

Supplemental Report:

No

CPT Code Information:

- 80358
- G0480 (if appropriate)

Reference Values:

Negative

Cutoff concentrations:

METHADONE BY GC-MS

<100 ng/mL

2-ETHYLIDENE-1,5-DIMETHYL-3,3-DIPHENYLPYRROLIDINE GC-MS

<100 ng/mL