

Reporting Title: Phenylalanine and Tyrosine, P
Performing Location: Rochester

- Necessary Information:
- 1. Patient's age is required.
 - 2. Include family history, clinical condition (asymptomatic or acute episode), diet, and drug therapy information.

- Specimen Requirements:
- Patient Preparation: Fasting (4 hours or more for infants)
- Collection Container/Tube:
- Preferred: Green top (sodium heparin)
- Acceptable: Green top (lithium heparin), Lavender top (EDTA)
- Submission Container/Tube: Plastic vial
- Specimen Volume: 0.5 mL
- Collection Instructions:
- 1. Centrifuge specimen and aliquot plasma into plastic vial.
 - 2. Send plasma frozen.

Forms:

[If not ordering electronically, complete, print, and send a Biochemical Genetics Test Request](#) (T798) with the specimen.

Specimen Type	Temperature	Time	Special Container
Plasma	Frozen (preferred)	14 days	
	Refrigerated	14 days	

Result Codes:

Result ID	Reporting Name	Type	Unit	LOINC®
8380	Phenylalanine, P	Numeric	nmol/mL	14875-9
8627	Tyrosine, P	Numeric	nmol/mL	20660-7

LOINC® and CPT codes are provided by the performing laboratory.

Supplemental Report:

No

- CPT Code Information:
- 84030-Phenylalanine
 - 84510-Tyrosine
 - 82542 (if appropriate for government payers)

Reference Values:

PHENYLALANINE

Premature: 98-213 nmol/mL

0-31 days: 38-137 nmol/mL
1-24 months: 31-75 nmol/mL
2-18 years: 26-91 nmol/mL
> or =19 years: 35-85 nmol/mL

Conversion Formulas:
Result in mg/dL x 60.5=result in nmol/mL
Result in nmol/mL x 0.0165=result in mg/dL

TYROSINE
Premature: 147-420 nmol/mL
0-31 days: 55-147 nmol/mL
1-24 months: 22-108 nmol/mL
2-18 years: 24-115 nmol/mL
> or =19 years: 34-112 nmol/mL

Conversion Formulas:
Result in mg/dL x 55.2=result in nmol/mL
Result in nmol/mL x 0.0181=result in mg/dL