

**Reporting Title:** Drug Screen, Prescription/OTC, S  
**Performing Location:** Rochester

**Ordering Guidance:**  
This test is not performed using chain of custody. For chain-of-custody testing, order DSSX / Drug Screen, Prescription/Over the Counter, Chain of Custody, Serum.

**Specimen Requirements:**  
**Supplies:** Sarstedt Aliquot Tube, 5 mL (T914)  
**Collection Container/Tube:** Red top (serum gel/SST tubes are **not** acceptable)  
**Submission Container/Tube:** Plastic vial  
**Specimen Volume:** 2.75 mL  
**Collection Instructions:** Centrifuge and aliquot serum into plastic vial within 2 hours of collection.

**Forms:**  
If not ordering electronically, complete, print, and send a [Therapeutics Test Request](#) (T831) with the specimen.

Specimen Type	Temperature	Time	Special Container
Serum Red	Refrigerated (preferred)	14 days	
	Frozen	14 days	
	Ambient	3 hours	

Ask at Order Entry (AOE) Questions:

Test ID	Question ID	Description	Type	Reportable
DSS	SU_DR	Suspect Drug	Plain Text	Yes

Result Codes:

Result ID	Reporting Name	Type	Unit	LOINC®
31168	Chain of Custody	Alphanumeric		77202-0
31072	Drugs detected:	Alphanumeric		20785-2

LOINC® and CPT codes are provided by the performing laboratory.

**Supplemental Report:**  
No

**CPT Code Information:**  
80307

Reference Values:

Drugs detected are presumptive. Additional testing may be required to confirm the presence of any drugs detected.