

T4 (Thyroxine), Total Only, Serum

**Reporting Title:** T4 (Thyroxine), Total Only, S **Performing Location:** Rochester

## **Ordering Guidance:**

This test cannot be used in patients receiving treatment with lipid-lowering agents containing dextrothyroxine unless therapy is discontinued for 4 to 6 weeks to allow the physiological state to become reestablished prior to testing.

## Specimen Requirements:

Collection Container/Tube: Preferred: Serum gel Acceptable: Red top Submission Container/Tube: Plastic vial Specimen Volume: 1 mL Collection Instructions: Centrifuge and aliquot serum into a plastic vial.

#### Forms:

Specimen Type	Temperature	Time	Special Container
Serum	Refrigerated (preferred)	7 days	
	Frozen	30 days	

### Result Codes:

Result ID	Reporting Name	Туре	Unit	LOINC®
T4	T4 (Thyroxine), Total Only, S	Numeric	mcg/dL	83119-8

 $\ensuremath{\mathsf{LOINC}}\xspace^{\ensuremath{\$}}$  and CPT codes are provided by the performing laboratory.

### Supplemental Report:

No

**CPT Code Information:** 84436

### **Reference Values:**

Pediatric 0-5 days: 5.0-18.5 mcg/dL 6 days-2 months: 5.4-17.0 mcg/dL 3-11 months: 5.7-16.0 mcg/dL 1-5 years: 6.0-14.7 mcg/dL 6-10 years: 6.0-13.8 mcg/dL 11-19 years: 5.9-13.2 mcg/dL



# **Test Definition: T4**

T4 (Thyroxine), Total Only, Serum

Adult (> or =20 years): 4.5-11.7 mcg/dL

For SI unit Reference Values, see <u>www.mayocliniclabs.com/order-tests/si-unit-conversion.html</u>