

**Reporting Title:** HBs Ab Monitor, Post-transplant, S  
**Performing Location:** Rochester

**Ordering Guidance:**

**Necessary Information:**  
Date of collection is required.

**Specimen Requirements:**

**Patient Preparation:** For 24 hours before specimen collection, patient **should not** take multivitamins or dietary supplements (eg, hair, skin, and nail supplements) containing biotin (vitamin B7).

**Supplies:** Sarstedt Aliquot Tube, 5 mL (T914)

**Collection Container/Tube:** Serum gel

**Submission Container/Tube:** Plastic vial

**Specimen Volume:** 0.7 mL

**Collection Instructions:**

1. Centrifuge blood collection tube per manufacturer's instructions (eg, centrifuge and aliquot within 2 hours of collection for BD Vacutainer tubes).
2. Aliquot serum into plastic vial.

**Forms:**

If not ordering electronically, complete, print, and send 1 of the following:

[-Gastroenterology and Hepatology Test Request](#) (T728)

[-Infectious Disease Serology Test Request](#) (T916)

Specimen Type	Temperature	Time	Special Container
Serum SST	Frozen (preferred)	90 days	
	Refrigerated	6 days	
	Ambient	72 hours	

**Result Codes:**

Result ID	Reporting Name	Type	Unit	LOINC®
HBABT	HBs Ab Monitor, Post-transplant, S	Alphanumeric	mIU/mL	5193-8

LOINC® and CPT codes are provided by the performing laboratory.

**Supplemental Report:**

No

**CPT Code Information:**

86317

Test Definition: HBABT

Hepatitis B Virus Surface Antibody Monitor,  
Post-Transplant, Serum

Reference Values:

Not applicable