

Reporting Title: Voriconazole, S

Performing Location: Rochester

Specimen Requirements:

Collection Container/Tube: Red top (serum gel/SST not acceptable)

Submission Container/Tube: Plastic vial

Specimen Volume: 2 mL

Collection Instructions: Centrifuge and aliquot serum into plastic vial within 2 hours of collection.

Forms:

| Specimen Type | Temperature              | Time    | Special Container |
|---------------|--------------------------|---------|-------------------|
| Serum Red     | Refrigerated (preferred) | 28 days |                   |
|               | Ambient                  | 28 days |                   |
|               | Frozen                   | 28 days |                   |

Result Codes:

| Result ID | Reporting Name  | Type    | Unit   | LOINC®  |
|-----------|-----------------|---------|--------|---------|
| 88698     | Voriconazole, S | Numeric | mcg/mL | 38370-3 |

LOINC® and CPT codes are provided by the performing laboratory.

Supplemental Report:

No

CPT Code Information:

80285

Reference Values:

1.0-5.5 mcg/mL

Trough level (ie, immediately before next dose) monitoring is recommended.