

Test Definition: C4U

C4 Acylcarnitine, Quantitative, Random, Urine

Reporting Title: C4 Acylcarnitine, QN, U

Performing Location: Rochester

Ordering Guidance:

This second-tier test is used specifically to evaluate a newborn screening elevation of iso-/butyrylcarnitine and **must not** be ordered with either C5OHU / C5-OH Acylcarnitine, Quantitative, Random, Urine or C5DCU / C5-DC Acylcarnitine, Quantitative, Random, Urine.

For general screening for metabolic disorders, see OAU / Organic Acids Screen, Random, Urine; ACRN / Acylcarnitines, Quantitative, Plasma; and AAQP / Amino Acids, Quantitative, Plasma.

Necessary Information:

Include patient's age, family history, clinical condition (asymptomatic or acute episode), diet, and drug therapy information.

Specimen Requirements:

Patient Preparation: If clinically feasible, discontinue L-carnitine supplementation at least 72 hours before specimen

collection.

Supplies: Urine Tubes, 10 mL (T068)

Collection Container/Tube: Clean, plastic urine collection container

Submission Container/Tube: Plastic, 10 mL urine tube

Specimen Volume: 5 mL **Collection Instructions:**

1. Collect a random urine specimen.

2. Freeze specimen immediately.

Forms:

If not ordering electronically, complete, print, and send a Biochemical Genetics Test Request (T798) with the specimen.

Specimen Type	Temperature	Time	Special Container
Urine	Frozen (preferred)	7 days	
	Refrigerated	24 hours	

Result Codes:

Result ID	Reporting Name	Туре	Unit	LOINC®
88829	C4 Acylcarnitine, QN, U	Numeric	mmol/mol Cr	53111-1
28075	C4 Interpretation	Alphanumeric		59462-2
34468	Reviewed By	Alphanumeric		18771-6

LOINC® and CPT codes are provided by the performing laboratory.

Supplemental Report:

No



Test Definition: C4U

C4 Acylcarnitine, Quantitative, Random, Urine

_	_	_	_				•								
-	u		T -	\sim	М.	$\hat{}$	\sim	_		m	•	+:	\sim	и	
	Ρ.			LJL		_			,, ,		_		.,		
•	•	•	•	•	ж.	•		•			•	•	•		

82017

Reference Values:

<3.00 millimoles/mole creatinine