

Reporting Title: Cyclobenzaprine (Flexeril)
Performing Location: Medtox Laboratories, Inc.

Specimen Requirements:
Submit only 1 of the following specimens:

Plasma
Specimen Type: Plasma
Container/Tube: Green Top
Specimen Volume: 2 mL
Collection Instructions: Draw blood in a green-top (sodium heparin) tube, **plasma gel tube is not acceptable**. Spin down and send 2 mL of sodium heparin plasma refrigerated in a plastic vial.

Serum
Specimen Type: Serum
Container/Tube: Red
Specimen Volume: 2 mL
Collection Instructions: Draw blood in a plain, red-top tube, **serum gel tube is not acceptable**. Spin down and send 2 mL of serum refrigerated in a plastic vial.

Specimen Type	Temperature	Time	Special Container
Varies	Refrigerated (preferred)	7 days	
	Frozen	30 days	
	Ambient	72 hours	

Result Codes:

Result ID	Reporting Name	Type	Unit	LOINC®
Z1130	Cyclobenzaprine (Flexeril)	Alphanumeric		6797-5

LOINC® and CPT codes are provided by the performing laboratory.

Supplemental Report:
No

CPT Code Information:
80369

Reference Values:
Reference Range: 10 - 30 ng/mL