

Reporting Title: Neurotensin
Performing Location: Inter Science Institute

Specimen Requirements:
Patient preparation: Patient should be fasting 10-12 hours prior to collection. Patient should not be on any antacid medication or medications that affect gastroentero-intestinal function, if possible, for at least 48 hours prior to collection.

Specimen Type: Plasma
Container/Tube: Z tube
Specimen Volume: 3 mL

Collection Instructions: Draw 10 mL of blood in special Z-tube (MCL T701). Separate plasma form cells immediately after draw and send 3 mL of plasma frozen in plastic vial.

Specimen Type	Temperature	Time	Special Container
Plasma	Frozen	30 days	

Result Codes:

Result ID	Reporting Name	Type	Unit	LOINC®
Z0923	Neurotensin	Numeric	pg/mL	2655-9

LOINC® and CPT codes are provided by the performing laboratory.

Supplemental Report:
No

CPT Code Information:
83519

Reference Values:
50 - 100 pg/mL