

Test Definition: FLIPR

Lipase, Random Urine

Reporting Title: Lipase, Urine

Performing Location: Quest Diagnostics Nichols Institute

Specimen Requirements:

Container/Tube: Plastic urine container

Specimen Volume: 2 mL Collection Instructions:

1. Collect random urine without preservative.

2. Ship 2 mL urine in a sterile screw capped plastic container ship ambient.

Specimen Type	Temperature	Time	Special Container	
Urine	Ambient (preferred)	5 days		
	Frozen	21 days		
	Refrigerated	5 days		

Result Codes:

Result ID	Reporting Name	Туре	Unit	LOINC®
Z1559	Lipase, Urine	Alphanumeric		10888-6

LOINC® and CPT codes are provided by the performing laboratory.

Supplemental Report:

No

CPT Code Information:

83690

Reference Values:

Adult: < or = 4 U/L