

# **Test Definition: FAMP**

Amphetamine, Serum or Plasma

Reporting Title: Amphetamine, Serum or Plasma

Performing Location: Medtox Laboratories, Inc.

## **Specimen Requirements:**

Submit only 1 of the following specimens:

#### Plasma

Draw blood in a green-top (sodium heparin) tube(s), plasma gel tube is not acceptable. Spin down and send 5 mL sodium heparin plasma refrigerated in a plastic vial.

#### Serum

Draw blood in a plain red-top tube(s), serum gel tube is not acceptable. Spin down and send 5 mL of serum refrigerated in a plastic vial.

Specimen Type	Temperature	Time	Special Container
Varies	Refrigerated (preferred)	14 days	
	Frozen	180 days	

### **Result Codes:**

Result ID	Reporting Name	Туре	Unit	LOINC®
Z1206	Amphetamine	Alphanumeric		42576-9

LOINC® and CPT codes are provided by the performing laboratory.

# **Supplemental Report:**

No

### **CPT Code Information:**

80324

# **Reference Values:**

Reference Range: 10 - 100 ng/mL