

**Reporting Title:** Pepsinogen I  
**Performing Location:** Inter Science Institute

**Specimen Requirements:**

**Patient preparation:**  
Patient should be fasting 10-12 hours prior to collection of specimen. Antacids or other medications affecting stomach acidity or gastrointestinal motility should be discontinued, if possible, for at least 48 hours prior to collection.

**Specimen Type:** Serum

**Container/Tube:** Red top or SST

**Specimen Volume:** 3 mL

**Collection Instructions:** Draw blood in a plain, red-top tube(s), serum-gel tube(s) is acceptable. Separate immediately and send 3 mL of serum frozen in a plastic vial.

Specimen Type	Temperature	Time	Special Container
Serum	Frozen	14 days	

**Result Codes:**

Result ID	Reporting Name	Type	Unit	LOINC®
Z0941	Pepsinogen I	Alphanumeric	ng/mL	2736-7

LOINC® and CPT codes are provided by the performing laboratory.

**Supplemental Report:**

No

**CPT Code Information:**

83520

**Reference Values:**

24 - 214 ng/mL