

Reporting Title: GlycoMark
Performing Location: Esoterix Endocrinology

Specimen Requirements:
Submit only 1 of the following specimens:

Serum
Draw blood in a plain red-top tube(s), serum gel tube is acceptable. Spin down and send 1 mL of serum ambient in a plastic vial.

Plasma
Draw blood in a lavender-top (EDTA) tube(s). Spin down and send 1 mL of EDTA plasma ambient in a plastic vial.

Specimen Type	Temperature	Time	Special Container
Varies	Ambient (preferred)	7 days	
	Frozen	14 days	
	Refrigerated	14 days	

Result Codes:

Result ID	Reporting Name	Type	Unit	LOINC®
Z0275	GlycoMark	Alphanumeric		53835-5

LOINC® and CPT codes are provided by the performing laboratory.

Supplemental Report:
No

CPT Code Information:
84378

Reference Values:

Age	Range
<18 y	Not Established
Adult Males	10.7-32.0
Adult Females	6.8-29.3

Glycemic control goal for diabetic patients: >10

GlycoMark is intended for use with managing glycemic control in diabetic patients. A low result corresponds to high glucose peaks.
1, 5-AG blood levels can be affected by clinical conditions or medications.