

Reporting Title: Flunitrazepam (Rohypnol)

Performing Location: Medtox Laboratories, Inc.

Specimen Requirements:

Submit only 1 of the following specimens:

Plasma

Draw blood in a green-top (sodium heparin) tube(s), **plasma gel tube is not acceptable**. Spin down and send 5 mL sodium heparin plasma refrigerated in a plastic vial.

Serum

Draw blood in a plain red-top tube(s), **serum gel tube is not acceptable**. Spin down and send 5 mL of serum refrigerated in a plastic vial.

Specimen Type	Temperature	Time	Special Container
Varies	Refrigerated (preferred)	14 days	
	Frozen	365 days	
	Ambient	48 hours	

Result Codes:

Result ID	Reporting Name	Type	Unit	LOINC®
Z1075	Flunitrazepam	Alphanumeric		73828-6
Z1074	7-Aminoflunitrazepam	Alphanumeric		73829-4

LOINC® and CPT codes are provided by the performing laboratory.

Supplemental Report:

No

CPT Code Information:

80346

Reference Values:

Units: Flunitrazepam ng/mL
 7-Aminoflunitrazepam ng/mL

Peak plasma Flunitrazepam concentrations in patients receiving chronic, recommended dosages: 10-20 ng/mL.

Note: Flunitrazepam is not legally marketed in the United States.