

Overview

Method Name
ImmunoCAP FEIA

NY State Available
Yes

Specimen

Specimen Type
Serum

Specimen Required
Draw blood in a plain, red-top tube(s), serum gel tube(s) is acceptable. Spin down and send 0.5 mL of serum refrigerated in a plastic vial.

Specimen Minimum Volume
0.5 mL

Reject Due To

Hemolysis:	NA
Thawing:	Warm OK; Cold OK
Lipemia:	Mild OK; Gross reject
Icterus:	NA
Other:	NA

Specimen Stability Information

Specimen Type	Temperature	Time	Special Container
Serum	Refrigerated (preferred)	28 days	
	Frozen	365 days	
	Ambient	28 days	

Clinical & Interpretive

Reference Values

<0.35 kU/L

Interpretation

Class	IgE (kU/L)	Comment
0	<0.10	Negative
0/1	0.10-0.34	Equivocal
1	0.35-0.69	Low Positive
2	0.70-3.4	Moderate Positive
3	3.5-17.4	High Positive
4	17.5-49.9	Very High Positive
5	50.0-99.9	Very High Positive
6	>100	Very High Positive

Performance

Method Description

The test method is the Phadia ImmunoCAP allergen-specific IgE system.

PDF Report

No

Day(s) Performed

Monday through Friday

Report Available

1 to 6 days

Performing Laboratory Location

Eurofins Viracor

Fees & Codes

Fees

- Authorized users can sign in to [Test Prices](#) for detailed fee information.
- Clients without access to Test Prices can contact [Customer Service](#) 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact [Customer Service](#).

Test Classification

[This test was developed and its performance characteristics determined by Viracor Eurofins. It has not been cleared or](#)

Test Definition: FWALP

Wall Eyed Pike (Sander vitreus)(Stizostedium vitreum) IgE

[approved by the U.S. Food and Drug Administration.](#)

CPT Code Information

86003

LOINC® Information

Test ID	Test Order Name	Order LOINC® Value
FWALP	Wall-Eyed Pike IgE	7761-0

Result ID	Test Result Name	Result LOINC® Value
Z3120	Wall-Eyed Pike IgE	7761-0
Z3121	Class	102796-0