

Overview

Method Name
Enzymatic/Spectrophotometric

NY State Available
Yes

Specimen

Specimen Type
Fecal

Specimen Required
Collect 1 gm random stool in sterile leak proof container, ship refrigerate.

Note:
Dietary restrictions: Patients receiving pancreatic enzymes should discontinue taking the enzymes at least 5 days before the collection of the stool sample.

Specimen Minimum Volume
0.5 gram

Reject Due To

Hemolysis	NA
Lipemia	NA
Icterus	NA
Other	NA

Specimen Stability Information

Specimen Type	Temperature	Time	Special Container
Fecal	Refrigerated (preferred)	14 days	
	Frozen	30 days	

Clinical & Interpretive

Reference Values
2.3 – 51.4 U/g

Performance

PDF Report

No

Day(s) Performed

Wednesday

Report Available

2 to 18 days

Performing Laboratory Location

Quest Diagnostics Nichols Institute

Fees & Codes

Fees

- Authorized users can sign in to [Test Prices](#) for detailed fee information.
- Clients without access to Test Prices can contact [Customer Service](#) 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact [Customer Service](#).

CPT Code Information

84311

LOINC® Information

Test ID	Test Order Name	Order LOINC® Value
FCHYS	Chymotrypsin, Stool	25375-7

Result ID	Test Result Name	Result LOINC® Value
FCHYS	Chymotrypsin, Stool	25375-7