

**Overview****Method Name**

ImmunoCAP FEIA

**NY State Available**

Yes

**Specimen****Specimen Type**

Serum

**Specimen Required**

Draw blood in a plain red-top tube(s), serum gel tube(s) is acceptable. Spin down and send 0.5 mL of serum refrigerated in a plastic vial.

**Specimen Minimum Volume**

0.5 mL

**Reject Due To**

Hemolysis	Mild OK; Gross reject
Lipemia	Mild OK; Gross reject
Icterus	Mild OK; Gross reject
Other	NA

**Specimen Stability Information**

Specimen Type	Temperature	Time	Special Container
Serum	Refrigerated (preferred)	28 days	
	Frozen	365 days	
	Ambient	28 days	

**Clinical and Interpretive****Reference Values**

&lt;0.35 kU/L

**Interpretation**

Class IgE (kU/L) Comment

0 &lt;0.10 Negative

0/1 0.10 – 0.34 Equivocal/Borderline

1 0.35 – 0.69 Low Positive

2 0.70 – 3.49 Moderate Positive

3 3.50 – 17.49 High Positive

4 17.50 – 49.99 Very High Positive

5 50.00 – 99.99 Very High Positive

6 >99.99 Very High Positive

### Performance

#### PDF Report

No

#### Day(s) Performed

Monday through Friday

#### Report Available

1 to 6 days

#### Performing Laboratory Location

Eurofins Viracor

### Fees and Codes

#### Fees

- Authorized users can sign in to [Test Prices](#) for detailed fee information.
- Clients without access to Test Prices can contact [Customer Service](#) 24 hours a day, seven days a week.
- Prospective clients should contact their Regional Manager. For assistance, contact [Customer Service](#).

#### CPT Code Information

86003

#### LOINC® Information

Test ID	Test Order Name	Order LOINC Value
FPRTF	Parrot Australian Feathers IgE	6030-1

Result ID	Test Result Name	Result LOINC Value
Z4052	Parrot Australian Feathers IgE	6030-1
Z4053	CLASS	15592-9