

Overview

Method Name

Hemolysis

NY State Available

Yes

Specimen

Specimen Type

Serum Red

Specimen Required

Specimen Type: Serum

Container/Tube: Red

Specimen Volume: 3 mL

Draw blood in a plain red-top tube. Maintain specimen at 37 degrees C until serum is separated from cells. Spin down and send 3 mL of serum refrigerated in a plastic vial.

Specimen Minimum Volume

2 mL

Reject Due To

|           |                       |
|-----------|-----------------------|
| Hemolysis | Mild OK; Gross reject |
| Lipemia   | Mild OK; Gross reject |
| Icterus   | Mild OK; Gross reject |
| Other     | Separator or gel tube |

Specimen Stability Information

| Specimen Type | Temperature              | Time     | Special Container |
|---------------|--------------------------|----------|-------------------|
| Serum Red     | Refrigerated (preferred) | 7 days   |                   |
|               | Ambient                  | 72 hours |                   |

Clinical & Interpretive

Reference Values

Negative

Performance

PDF Report

No

Day(s) Performed

Monday through Friday

Report Available

1 to 7 days

Performing Laboratory Location

ARUP Laboratories

Fees & Codes

Fees

- Authorized users can sign in to [Test Prices](#) for detailed fee information.
- Clients without access to Test Prices can contact [Customer Service](#) 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact [Customer Service](#).

CPT Code Information

86940, 86941

LOINC® Information

| Test ID | Test Order Name    | Order LOINC® Value |
|---------|--------------------|--------------------|
| FDLS    | Donath Landsteiner | 18287-3            |

| Result ID | Test Result Name   | Result LOINC® Value |
|-----------|--------------------|---------------------|
| FDLS      | Donath Landsteiner | 18287-3             |