Test Definition: FGHSP

GHB Screen, S/P

Overview

Reflex Tests

Test ID	Reporting Name	Available Separately	Always Performed
FGHA	Gamma-Hydroxybutyric Acid, CF, SP	No	No

Testing Algorithm

If the GHB Screen S/P is non-negative, then the Gamma-Hydroxybutyric Acid, CF, SP (FGHA) will be performed at an additional charge.

Method Name

Liquid Chromatography with Tandem Mass Spectrometry (LC/MS/MS)

Gas Chromatography/Mass Spectrometry (GC/MS) (if appropriate)

NY State Available

Yes

Specimen

Specimen Type

Varies

Specimen Required

Submit only 1 of the following specimens:

Serum

Draw blood in a plain red-top tube(s), **serum gel tube is not acceptable**. Spin down and send 5 mL of serum refrigerated in a plastic vial.

Plasma

Draw blood in a green-top (sodium heparin) tube(s), **plasma gel tube is not acceptable**. Spin down and send 5 mL of sodium heparin plasma refrigerated in a plastic vial.

Specimen Minimum Volume

1.2 mL

Reject Due To

Hemolysis	NA
Lipemia	NA
Icterus	NA



Test Definition: FGHSP

GHB Screen, S/P

Other	NA

Specimen Stability Information

Specimen Type	Temperature	Time	Special Container
Varies	Refrigerated (preferred)	7 days	
	Frozen	180 days	
	Ambient	72 hours	

Clinical and Interpretive

Reference Values

Reference Range: Negative

Screening threshold: 5.0 ug/mL

Performance

PDF Report

No

Day(s) and Time(s) Test Performed

Monday through Sunday

Analytic Time

5 days

Maximum Laboratory Time

7 - 9 days

Performing Laboratory Location

Medtox Laboratories, Inc.

Fees and Codes

Fees

- Authorized users can sign in to <u>Test Prices</u> for detailed fee information.
- Clients without access to Test Prices can contact <u>Customer Service</u> 24 hours a day, seven days a week.
- Prospective clients should contact their Regional Manager. For assistance, contact Customer Service.

CPT Code Information

80307

LOINC® Information



Test Definition: FGHSP

GHB Screen, S/P

Test ID	Test Order Name	Order LOINC Value
FGHSP	GHB Screen, S/P	46085-7

Result ID	Test Result Name	Result LOINC Value
Z4223	GHB	46085-7