

Overview

Method Name

Only orderable as part of a profile. For more information see:

- MSP3 / Multiple Sclerosis (MS) Profile, Serum and Spinal Fluid
- OLIG / Oligoclonal Banding, Serum and Spinal Fluid

Collection Only

NY State Available

No

Specimen

Specimen Type

Serum

Specimen Required

Only orderable as part of a profile. For more information see:

- MSP3 / Multiple Sclerosis (MS) Profile, Serum and Spinal Fluid
- OLIG / Oligoclonal Banding, Serum and Spinal Fluid

Specimen Type: Serum

Collection Container/Tube:

Preferred: Serum gel

Acceptable: Red top

Submission Container/Tube: Plastic vial

Specimen Volume: 1 mL

Collection Instructions:

1. Centrifuge and aliquot serum in plastic vial within 2 hours of collection.
2. Label specimen as serum.

Specimen Minimum Volume

0.5 mL

Reject Due To

Gross hemolysis	OK
Gross lipemia	OK

Gross icterus	OK
---------------	----

Specimen Stability Information

Specimen Type	Temperature	Time	Special Container
Serum	Frozen (preferred)	14 days	
	Ambient	14 days	
	Refrigerated	14 days	

Clinical & Interpretive

Performance

PDF Report
No

Specimen Retention Time
14 days

Performing Laboratory Location
Rochester

Fees & Codes

Fees

- Authorized users can sign in to [Test Prices](#) for detailed fee information.
- Clients without access to Test Prices can contact [Customer Service](#) 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact [Customer Service](#).

Test Classification
Not Applicable

LOINC® Information

Test ID	Test Order Name	Order LOINC® Value
XSRM	Additional sample for reflex OLIGS	No LOINC Needed

Result ID	Test Result Name	Result LOINC® Value
XSRM	Additional sample for Reflex OLIGS	No LOINC Needed