

Overview**Method Name**

Liquid Chromatography/Tandem Mass Spectrometry (LC/MS/MS)

NY State Available

Yes

Specimen**Specimen Type**

Varies

Specimen Required

Submit only 1 of the following specimens:

Serum

Draw blood in a plain, red-top tube(s), serum gel tube is not acceptable. Spin down and send 2 mL of serum refrigerated in a plastic vial.

Plasma

Draw blood in a green-top sodium heparin tube(s), plasma gel tube is not acceptable. Spin down and send 2 mL of plasma refrigerated in a plastic vial.

Specimen Minimum Volume

1.0 mL

Reject Due To

Hemolysis	NA
Lipemia	NA
Icterus	NA
Other	Serum gel, Plasma gel

Specimen Stability Information

Specimen Type	Temperature	Time	Special Container
Varies	Refrigerated (preferred)	14 days	
	Frozen	180 days	
	Ambient	72 hours	

Clinical and Interpretive

Reference Values

5 - 25 ng/mL

Reporting Limit: 2.0 ng/mL

Performance**PDF Report**

No

Day(s) and Time(s) Test Performed

Monday through Sunday

Analytic Time

5 days

Maximum Laboratory Time

7 - 9 days

Performing Laboratory Location

Medtox Laboratories, Inc.

Fees and Codes**Fees**

- Authorized users can sign in to [Test Prices](#) for detailed fee information.
- Clients without access to Test Prices can contact [Customer Service](#) 24 hours a day, seven days a week.
- Prospective clients should contact their Regional Manager. For assistance, contact [Customer Service](#).

Test Classification

In-house validated method

CPT Code Information

80346

LOINC® Information

Test ID	Test Order Name	Order LOINC Value
FALPX	Alprazolam (Xanax)	59611-4

Result ID	Test Result Name	Result LOINC Value
FALPX	Alprazolam (Xanax)	59611-4