

Overview

Special Instructions

- [National Jewish Complement request form](#)

Method Name

Radioimmunoassay (RIA)

NY State Available

Yes

Specimen

Specimen Type

Plasma EDTA

Specimen Required

Container/Tube: Lavender top tube

Preferred: Lavender top tube

Acceptable: Lavender top tube

Specimen Volume: 2 mL

Collection Instructions: Draw blood in a lavender-top (EDTA) tube(s), **plasma gel tube(s) is not acceptable**. Mix well, centrifuge at ambient temperature within one half hour of draw and freeze immediately. Send 1 mL of EDTA plasma frozen in plastic vial.

Complete and submit with specimen:

[National Jewish Complement request form](#)

Specimen Minimum Volume

0.5 mL

Reject Due To

Hemolysis	NA
Lipemia	NA
Icterus	NA
Other	NA

Specimen Stability Information

Specimen Type	Temperature	Time	Special Container
---------------	-------------	------	-------------------

Plasma EDTA	Frozen	365 days	
-------------	--------	----------	--

Clinical & Interpretive

Reference Values

0-780 ng/mL

Performance

PDF Report

No

Day(s) Performed

First and third Thursday of the month

Report Available

21 to 44 days

Performing Laboratory Location

National Jewish Health

Fees & Codes

Fees

- Authorized users can sign in to [Test Prices](#) for detailed fee information.
- Clients without access to Test Prices can contact [Customer Service](#) 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact [Customer Service](#).

Test Classification

The performance characteristics of this test have been validated by National Jewish Health. It has not been cleared or approved by the U.S. Food and Drug Administration. The results are not intended to be used as the sole means for clinical diagnosis or patient management decisions. This laboratory is certified under the Clinical Laboratory Improvement Amendments of 1988 (CLIA-88) as qualified to perform high complexity clinical laboratory testing.

CPT Code Information

86160

LOINC® Information

Test ID	Test Order Name	Order LOINC® Value
FC3AR	C3a Level	4488-3

Result ID	Test Result Name	Result LOINC® Value
FC3AR	C3a Level	4488-3