



# Test Definition: FASQN

Antidepressant Drug Screen, Ur, Quantitative

## Overview

### Method Name

Immunoassay (IA); Liquid Chromatography/Tandem Mass Spectrometry (LC/MS/MS)

### NY State Available

Yes

## Specimen

### Specimen Type

Urine

### Specimen Required

Collect 3 mL random urine without preservatives. Ship refrigerated in a plastic container.

### Specimen Minimum Volume

0.6 mL

### Reject Due To

Hemolysis	NA
Lipemia	NA
Icterus	NA
Other	NA

### Specimen Stability Information

Specimen Type	Temperature	Time	Special Container
Urine	Refrigerated (preferred)	14 days	
	Ambient	72 hours	
	Frozen	30 days	

## Clinical & Interpretive

### Reference Values

Antidepressant screen includes the analysis for: amitriptyline, clomipramine and desmethylclomipramine, cyclobenzaprine, desipramine, doxepin and desmethldoxepin, fluoxetine and norfluoxetine, imipramine, maprotiline, nortriptyline, paroxetine, protriptyline, sertraline and desmethylsertraline, trimipramine.

**Performance****PDF Report**

No

**Day(s) Performed**

Monday through Sunday

**Report Available**

3 to 9 days

**Performing Laboratory Location**

Medtox Laboratories, Inc.

**Fees & Codes****Fees**

- Authorized users can sign in to [Test Prices](#) for detailed fee information.
- Clients without access to Test Prices can contact [Customer Service](#) 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact [Customer Service](#).

**Test Classification**

This test was developed and its performance characteristics determined by LabCorp. It has not been cleared or approved by the Food and Drug Administration.

**CPT Code Information**

80333

80337

80369

**LOINC® Information**

Test ID	Test Order Name	Order LOINC® Value
FASQN	Antidepressant Screen, Quant, U	32136-4

Result ID	Test Result Name	Result LOINC® Value
Z3350	Antidepressant Screen	32136-4
Z3351	Amitriptyline	16225-5
Z3352	Nortriptyline	61427-1
Z3353	Imipramine	61418-0
Z3354	Desipramine	61412-3
Z3355	Sertraline	32112-5

Z3356	Desmethylsertraline	75233-7
Z3357	Fluoxetine	19471-2
Z3358	Norfluoxetine	19629-5
Z3359	Clomipramine	45144-3
Z3360	Desmethylclomipramine	58013-4
Z3361	Doxepin	61415-6
Z3362	Desmethyldoxepin	58014-2
Z3363	Maprotiline	3737-4
Z3364	Protriptyline	4001-4
Z3365	Trimipramine	4084-0
Z3366	Paroxetine	27109-8
Z3367	Cyclobenzaprine	14061-6