

Overview

Method Name
Enzyme Linked Immunosorbent Assay (ELISA)

NY State Available
Yes

Specimen

Specimen Type
Serum

Specimen Required
Specimen Type: Serum
Container/Tube: Red or SST
Specimen Volume:1 mL
Collection Instructions: Draw blood in a plain, red-top tube(s), **serum gel tube(s) is acceptable**. Spin down and send 1 mL frozen serum in a screw-capped vial.

Specimen Minimum Volume
0.2 mL

Reject Due To

Hemolysis	Mild reject; Gross reject
-----------	---------------------------

Specimen Stability Information

Specimen Type	Temperature	Time	Special Container
Serum	Frozen		

Clinical & Interpretive

Reference Values
Reference Range: 0.0 - 15.0 ng/mL

Note: Graph attached as a supplemental report in MayoAccess

Clinical Reference

Fehm T, Malmonis P, Weitz S, Teramoto Y, Katalinic A, Jager W. Influence of circulating c-erbB-2 serum protein on response to adjuvant chemotherapy in node-positive breast cancer patients. Breast Cancer Res Treat. 1997 Mar; 43(1):87-95.

Molina R, Jo J, Filella X, et al., c-erbB-2 oncoprotein in the sera and tissue of patients with breast cancer: Utility in prognosis. Anti-cancer Res. 1996 Jul-Aug; 16(4B) 2295-2300.

Performance

Method Description

Bayer Centaur/ACS methodology

PDF Report

Referral

Day(s) Performed

Tuesday, Friday

Report Available

7 to 9 days

Performing Laboratory Location

LabCorp Burlington

Fees & Codes

Fees

- Authorized users can sign in to [Test Prices](#) for detailed fee information.
- Clients without access to Test Prices can contact [Customer Service](#) 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact [Customer Service](#).

CPT Code Information

83950

LOINC® Information

Test ID	Test Order Name	Order LOINC® Value
FHER	HER-2/neu Quantitative ELISA	42914-2

Result ID	Test Result Name	Result LOINC® Value
Z0985	HER-2/neu Quantitative	42914-2