

# **Reporting Title:** Glucose, CSF **Performing Location:** Rochester

#### **Specimen Requirements:**

Container/Tube: Plastic vial Specimen Volume: 1 mL Collection Instructions: Centrifuge to remove any cellular material.

#### **Specimen Minimum Volume:**

0.25 mL

Specimen Type	Temperature	Time	Special Container
CSF	Frozen (preferred)	30 days	
	Refrigerated	7 days	

### Result Codes:

Result ID	Reporting Name	Туре	Unit	LOINC®
GLSF	Glucose, CSF	Numeric	mg/dL	2342-4

LOINC and CPT codes are provided by the performing laboratory.

#### Supplemental Report:

No

# **CPT Code Information:**

82945

## **Reference Values:**

Spinal fluid glucose concentration should be approximately 60% of the plasma/serum concentration and should be compared with concurrently measured plasma/serum glucose for adequate clinical interpretation.

For SI unit Reference Values, see https://www.mayocliniclabs.com/order-tests/si-unit-conversion.html