



Reporting Title: HBc Total Ab, w/Reflex, S

Performing Location: Rochester

Necessary Information:

Date of collection is required.

Specimen Requirements:

Patient Preparation: For 24 hours before specimen collection, patient should not take multivitamins or dietary supplements (eg, hair, skin, and nail supplements) containing biotin (vitamin B7).

Collection Container/Tube: Sarstedt Aliquot Tube, 5 mL (T914)

Submission Container/Tube: Plastic vial

Specimen Volume: 0.7 mL

Collection Instructions:

1. Centrifuge blood collection tube per manufacturer's instructions (eg, centrifuge and aliquot within 2 hours of collection for BD Vacutainer tubes).
2. Aliquot serum into a plastic vial.

Specimen Minimum Volume:

0.6 mL

Forms:

If not ordering electronically, complete, print, and send 1 of the following:

-Gastroenterology and Hepatology Test Request (T728)

-Infectious Disease Serology Test Request (T916)

Specimen Type	Temperature	Time	Special Container
Serum SST	Frozen (preferred)	90 days	
	Refrigerated	6 days	
	Ambient	72 hours	

Result Codes:

Result ID	Reporting Name	Type	Unit	LOINC®
CORAB	HBc Total Ab, w/Reflex, S	Alphanumeric		13952-7

LOINC and CPT codes are provided by the performing laboratory.



Supplemental Report:

No

CPT Code Information:

86704

86705 (if appropriate)

Reflex Tests:

Test ID	Reporting Name	CPT Units	CPT Code	Always Performed	Orderable Separately
HBIM	HBc IgM Ab, S			No	Yes

Result Codes for Reflex Tests:

Test ID	Result ID	Reporting Name	Type	Unit	LOINC®
HBIM	HBIM	HBc IgM Ab, S	Alphanumeric		24113-3

Reference Values:

Negative

Interpretation depends on clinical setting.

See Viral Hepatitis Serologic Profiles