

Tryptase, Autopsy, Serum

# **Reporting Title:** Tryptase, Autopsy **Performing Location:** Rochester

### **Specimen Requirements:**

Collection Container/Tube: Preferred: Serum gel Acceptable: Red top Submission Container/Tube: Plastic vial Specimen Volume: 0.5 mL Collection Instructions: Centrifuge and aliquot serum into plastic vial. Additional Information: Tryptase degenerates very quickly when left in the presence of red blood cells.

### **Specimen Minimum Volume:**

0.2 mL

Specimen Type	Temperature	Time	Special Container
Serum	Frozen (preferred)	14 days	
	Refrigerated	7 days	

## **Result Codes:**

Result ID	Reporting Name	Туре	Unit	LOINC®
TRYPA	Tryptase, Autopsy	Numeric	ng/mL	21582-2

LOINC and CPT codes are provided by the performing laboratory.

#### **Supplemental Report:**

No

### **CPT Code Information:**

83520



Tryptase, Autopsy, Serum

#### **Reference Values:**

No established reference values