
Reporting Title: Fragile X Syndrome, Mol. Analysis**Performing Location:** Rochester**Ordering Guidance:**

Due to the complexity of prenatal testing, consultation with the laboratory is required. To speak with a genetic counselor about this testing option, call 800-533-1710.

FMR1-methylation status cannot be assessed on chorionic villus specimens. Contact a molecular genetic counselor/consultant at 800-533-1710 to discuss the limitations of testing prior to sending a chorionic villus specimen for fragile X analysis.

Additional Testing Requirements:

All prenatal specimens must be accompanied by a maternal blood specimen. Order MATCC / Maternal Cell Contamination, Molecular Analysis, Varies on the maternal specimen. This must be a different order number than the prenatal specimen.

Shipping Instructions:

Specimen preferred to arrive within 96 hours of collection.

Specimen Requirements:

Patient Preparation: A previous bone marrow transplant from an allogenic donor will interfere with testing. For instructions for testing patients who have received a bone marrow transplant, call 800-533-1710.

Submit only 1 of the following specimens:

Specimen Type: Whole blood

Container/Tube:

Preferred: Lavender top (EDTA) or yellow top (ACD)

Acceptable: Any anticoagulant

Specimen Volume: 3 mL

Collection Instructions:

1. Invert several times to mix blood.

2. Send whole blood specimen in original tube. Do not aliquot.

Specimen Stability Information: Ambient (preferred) 4 days/Refrigerated

Prenatal Specimens

Due to its complexity, consultation with the laboratory is required for all prenatal testing; call 800-533-1710 to speak to a genetic counselor.

Specimen Type: Amniotic fluid

Container/Tube: Amniotic fluid container

Specimen Volume: 20 mL

Specimen Stability Information: Refrigerated (preferred)/Ambient

Additional information:

-
1. A separate culture charge will be assessed under CULAF / Culture for Genetic Testing, Amniotic Fluid. An additional 2 to 3 weeks is required to culture amniotic fluid before genetic testing can occur.
 2. All prenatal specimens must be accompanied by a maternal blood specimen; order MATCC / Maternal Cell Contamination, Molecular Analysis, Varies on the maternal specimen.

Specimen Type: Chorionic villi

Container/Tube: 15-mL tube containing 15 mL of transport media

Specimen Volume: 20 mg

Specimen Stability Information: Refrigerated

Additional Information:

1. A separate culture charge will be assessed under CULFB / Fibroblast Culture for Biochemical or Molecular Testing. An additional 2 to 3 weeks is required to culture chorionic villi before genetic testing can occur.
2. All prenatal specimens must be accompanied by a maternal blood specimen; order MATCC / Maternal Cell Contamination, Molecular Analysis, Varies on the maternal specimen.

Specimen Type: Confluent cultured cells

Container/Tube: T-25 flask

Specimen Volume: 2 Flasks

Collection Instructions: Submit confluent cultured cells from another laboratory.

Specimen Stability Information: Ambient (preferred)/Refrigerated

Additional Information: All prenatal specimens must be accompanied by a maternal blood specimen; order MATCC / Maternal Cell Contamination, Molecular Analysis, Varies on the maternal specimen.

Specimen Minimum Volume:

Blood: 0.5 mL

Amniotic fluid: 10 mL

Chorionic villi: 5 mg

Forms:

1. New York Clients-Informed consent is required. Document on the request form or electronic order that a copy is on file. The following documents are available:
 - Informed Consent for Genetic Testing (T576)
 - Informed Consent for Genetic Testing-Spanish (T826)
2. Molecular Genetics: Congenital Inherited Diseases Patient Information (T521)
3. If not ordering electronically, complete, print, and send a Neurology Specialty Testing Client Test Request (T732) with the specimen.

Specimen Type	Temperature	Time	Special Container
Varies	Varies		

Ask at Order Entry (AOE) Questions:

Test ID	Question ID	Description	Type	Reportable
CULFB	CG770	Reason for Referral	Plain Text	No
CULFB	CG899	Specimen	Plain Text	No
CULAF	CG767	Reason for Referral	Plain Text	No

Result Codes:

Result ID	Reporting Name	Type	Unit	LOINC®
52870	Result Summary	Alphanumeric		50397-9
52871	Result	Alphanumeric		81856-7
52872	Interpretation	Alphanumeric		69047-9
52873	Reason for Referral	Alphanumeric		42349-1
52874	Specimen	Alphanumeric		31208-2
52875	Source	Alphanumeric		31208-2
52876	Method	Alphanumeric		85069-3
52877	Released By	Alphanumeric		18771-6

LOINC and CPT codes are provided by the performing laboratory.

Supplemental Report:

No

CPT Code Information:

81243
88233 (if appropriate)
88240 (if appropriate)
88235 (if appropriate)
81265 (if appropriate)
81244 (if appropriate)

Reflex Tests:

Test ID	Reporting Name	CPT Units	CPT Code	Always Performed	Orderable Separately
CULFB	Fibroblast Culture for Genetic Test			No	Yes
CULAF	Amniotic Fluid Culture/Genetic Test			No	Yes
MATCC	Maternal Cell Contamination, B			No	Yes
_STR1	Comp Analysis using STR (Bill only)			No	No (Bill only)
_STR2	Add'l comp analysis w/STR (Bill Only)			No	No (Bill only)
FUFXS	Fragile X, Follow up Analysis			No	No

Result Codes for Reflex Tests:

Test ID	Result ID	Reporting Name	Type	Unit	LOINC®
CULFB	52327	Result Summary	Alphanumeric		50397-9
CULFB	52329	Interpretation	Alphanumeric		69965-2
CULFB	52328	Result	Alphanumeric		82939-0
CULFB	CG770	Reason for Referral	Alphanumeric		42349-1
CULFB	CG899	Specimen	Alphanumeric		31208-2
CULFB	52331	Source	Alphanumeric		31208-2
CULFB	52332	Method	Alphanumeric		85069-3
CULFB	54625	Additional Information	Alphanumeric		48767-8
CULFB	52333	Released By	Alphanumeric		18771-6
CULAF	52304	Result Summary	Alphanumeric		50397-9
CULAF	52306	Interpretation	Alphanumeric		69965-2
CULAF	52305	Result	Alphanumeric		82939-0
CULAF	CG767	Reason for Referral	Alphanumeric		42349-1
CULAF	52307	Specimen	Alphanumeric		31208-2
CULAF	52308	Source	Alphanumeric		31208-2
CULAF	52309	Method	Alphanumeric		85069-3
CULAF	54641	Additional Information	Alphanumeric		48767-8

Test ID	Result ID	Reporting Name	Type	Unit	LOINC®
CULAF	52310	Released By	Alphanumeric		18771-6
MATCC	53285	Result Summary	Alphanumeric		50397-9
MATCC	53286	Result	Alphanumeric		40704-9
MATCC	53287	Interpretation	Alphanumeric		69047-9
MATCC	53288	Reason for referral	Alphanumeric		42349-1
MATCC	53289	Specimen	Alphanumeric		31208-2
MATCC	53290	Source	Alphanumeric		31208-2
MATCC	55150	Method	Alphanumeric		85069-3
MATCC	53291	Released By	Alphanumeric		18771-6
FUFXS	52421	Comment	Alphanumeric		48767-8
FUFXS	52422	Specimen	Alphanumeric		31208-2
FUFXS	52423	Source	Alphanumeric		31208-2
FUFXS	52424	Released By	Alphanumeric		18771-6

Reference Values:

Normal alleles: 5-44 CGG repeats

Intermediate (grey zone) alleles: 45-54 CGG repeats

Premutation alleles: 55-200 CGG repeats

Full mutation alleles: >200 CGG repeats

An interpretive report will be provided.

Methylation status:

Unmethylated: < or =20%

Partially methylated: 21-69%

Fully methylated: > or =70%