

Valproic Acid, Total, Serum

# **Reporting Title:** Valproic Acid, Tot, S **Performing Location:** Rochester

### **Specimen Requirements:**

Collection Container/Tube: Preferred: Serum gel Acceptable: Red top Submission Container/Tube: Plastic vial Specimen Volume: 0.5 mL Collection Instructions: 1. Serum gel tubes should be centrifuged within 2 hours of collection. 2. Red-top tubes should be centrifuged and serum aliquoted into a plastic vial within 2 hours of collection.

### **Specimen Minimum Volume:**

0.25 mL

### Forms:

If not ordering electronically, complete, print, and send 1 of the following forms with the specimen: -Neurology Specialty Testing Client Test Request (T732) -Therapeutics Test Request (T831)

Specimen Type	Temperature	Time	Special Container
Serum	Refrigerated (preferred)	14 days	
	Frozen	28 days	
	Ambient	72 hours	

### **Result Codes:**

Result ID	Reporting Name	Туре	Unit	LOINC®
VALPA	Valproic Acid, Tot, S	Numeric	mcg/mL	4086-5

LOINC and CPT codes are provided by the performing laboratory.

## Supplemental Report:

No

### **CPT Code Information:**

80164



Valproic Acid, Total, Serum

### **Reference Values:**

Therapeutic: 50 (trough)-125 (peak) mcg/mL Critical value: > or =151 mcg/mL