

**Reporting Title:** Hepatitis A IgM Ab, S**Performing Location:** Rochester**Necessary Information:**

Date of collection is required.

**Specimen Requirements:**

Patient Preparation: For 24 hours before specimen collection, patient should not take multivitamins or dietary supplements (eg, hair, skin, and nail supplements) containing biotin (vitamin B7).

Collection Container/Tube: Serum gel

Submission Container/Tube: Plastic vial

Specimen Volume: 0.6 mL

Collection Instructions:

1. Centrifuge blood collection tube per manufacturer's instructions (eg, centrifuge and aliquot within 2 hours of collection for BD Vacutainer tubes).
2. Aliquot serum into plastic vial.

**Specimen Minimum Volume:**

0.6 mL

**Forms:**

If not ordering electronically, complete, print, and send 1 of the following:

-Gastroenterology and Hepatology Test Request (T728)

-Infectious Disease Serology Test Request (T916)

Specimen Type	Temperature	Time	Special Container
Serum SST	Frozen (preferred)	90 days	
	Refrigerated	6 days	
	Ambient	72 hours	

**Result Codes:**

Result ID	Reporting Name	Type	Unit	LOINC®
HAIGM	Hepatitis A IgM Ab, S	Alphanumeric		13950-1

LOINC and CPT codes are provided by the performing laboratory.

---

**Supplemental Report:**

No

**CPT Code Information:**

86709

**Reference Values:**

Negative

See Viral Hepatitis Serologic Profiles.