

Trofile DNA Co-Receptor Tropism Assay

Reporting Title: Trofile DNA Phenotypic Corecept Tro **Performing Location:** Monogram Biosciences, Inc

Specimen Requirements:

Specimen Type: Whole Blood Container/Tube: Lavender-top (EDTA) Specimen Volume: 4 mL Collection Instructions: Draw 4 mL blood in a lavender-top (EDTA) tube(s), (Do not centrifuge.) Freeze and ship frozen. To avoid delays in turnaround time when requesting multiple tests, please submit separate frozen specimens for each test requested. Note: Trofile DNA is recommended for patients with undetectable viral loads.

Specimen Minimum Volume:

3 mL

Specimen Type	Temperature	Time	Special Container
Whole Blood EDTA	Frozen	14 days	

Result Codes:

Result ID	Reporting Name	Туре	Unit	LOINC®
Z2283	Tropotype Result	Alphanumeric		53923-9
Z2284	Interpretation:	Alphanumeric		Not Provided

LOINC and CPT codes are provided by the performing laboratory.

Supplemental Report:

No

CPT Code Information:

87999



Test Definition: FFTRO

Trofile DNA Co-Receptor Tropism Assay

Reference Values: