

Reporting Title: Maple Red (Acer rubrum) IgE**Performing Location:** Eurofins Viracor**Specimen Requirements:**

Draw blood in a plain red-top tube(s), serum gel tube(s) is acceptable. Spin down and send 0.5 mL of serum refrigerated in a plastic vial.

Specimen Minimum Volume:

0.5 mL

Specimen Type	Temperature	Time	Special Container
Serum	Refrigerated (preferred)	28 days	
	Frozen	365 days	
	Ambient	28 days	

Result Codes:

Result ID	Reporting Name	Type	Unit	LOINC®
Z3094	Maple Red (Acer rubrum) IgE	Alphanumeric		7479-9
Z3095	Class	Alphanumeric		102427-2

LOINC and CPT codes are provided by the performing laboratory.

Supplemental Report:

No

CPT Code Information:

86003

Reference Values:

<0.35 kU/L